# ***Instructions for Use***

Instructions for Use - GB  
Gebrauchsanweisung - D  
Instructions d'utilisation - F  
Bruksanvisning - S  
Istruzioni per l'Uso - I  
Instrucções para a utilização - P  
Käyttöohje - FI  
Οδηγίες χρήσεως - GR  
Manual de instrucciones - E  
Brugsanvisning - DK  
Návod k použití - CZ  
Instrukcja obsługi - PL  
Naudotojo vadovas - LT  
Kasutusjuhend - EE  
Használati utasítás - HU  
Gebruiksaanwijzing - NL  
Instrucțiuni de utilizare - RO

## ***Impedance Audiometer AT235***

Valid from serial 639427  
80650305 – 07/2008

# ***Instructions for Use - GB***

## ***Impedance Audiometer AT235***

## Intended use

The AT235 is an automatic impedance audiometer suited for clinics performing screening, as well as diagnostic work. The test battery includes tympanometry, acoustic reflex testing, ipsilateral and contralateral reflex decay testing, Eustachian tube function test and AC audiometry.

The AT235 tympanometer is intended to be used by an audiologist, hearing healthcare professional, or trained technician in a quiet environment (tympanometry and reflexes) and extremely quiet environment (Audiometry). Careful handling of instrument whenever in contact with patient should be of high priority. Calm and stable positioning while testing is preferred for optimal accuracy. It is recommended that the instrument be operated within an ambient temperature range of 15-35 degree Celsius (59-95 degrees Fahrenheit).

## Precautions

**Notice** - Be sure to insert the probe tip in a way, which will assure an air tight fit without causing any harm to the patient. Using a proper and clean ear tip is mandatory.

**Notice** - We recommend using a new ear tip for each patient.

If the clinician rinses the ear tips they should be subjected to standard disinfecting procedure between patients. This includes physically cleaning the ear tip and use of a recognised disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriated level of cleanliness.

**Notice** - Be sure to use only stimulation intensities, which will be acceptable for the patient.

**Notice** - The transducers (headphones, bone conductor, etc.) supplied with the instrument are calibrated to this instrument - exchange of transducers require a re-calibration.

**Notice** - It is recommended that parts which are in direct contact with the patient (e.g. earphone cushions) are subjected to standard disinfecting procedure between patients. This includes physically cleaning and use of a recognised disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriated level of cleanliness.

**Notice** - Although the instrument fulfils the relevant EMC requirements precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g. from mobile phones etc. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears

**Notice** - Please note that the CE Marking is only legal if this Instruction is translated into the national language of the user no later than at the delivery to him, if the national legislation demands a text in the national language according to MDD article 4.4.

If this apparatus is connected to one or more other devices with medical CE marking, to make up a system or pack, the CE marking is only valid also for the combination if the supplier has issued a declaration stating that the requirements in the Medical Device Directive article 12 are fulfilled for the combination.

## Operation Instruction

### 1. Switch on the instrument

Make sure that the power supply is connected correctly to the mains and to the instrument. Now switch on the instrument on the connection panel of the instrument by pressing the power switch (33).

### 2. Clear memory of instrument

Before making a test on a new patient it is always recommended to clear the internal memory of the instrument by holding down "shift" (17) while pressing "new subject" (9).

### 3. The right eartip

Check the size of the patient's ear canal in order to find an appropriate size eartip (6 to 22 mm) and mount it onto the probe tip. It is recommended always to use the eartips supplied by Interacoustics.

### 4. Select ear

Now select the correct ear by pressing "Right" or "Left" (20 or 21) or by pressing the remote switch on the probe.

### 5. Testing the first ear

Take the probe with one hand and insert the probe tip with the eartip mounted into the patient's ear canal while, with the other hand, carefully pulling backwards / upwards the pinna to straighten out the ear canal during the insertion.

The probe will have an air tight fit when the indication light on the probe is green. When the status "Leaking" or "Blocked" is displayed in the upper right hand corner of the screen under the Status Bar, the indication light on the probe will be yellow. This indicates a bad fitting or a blocked ear tip.

When a correct fit has been obtained a tympanometric test is automatically performed on the first ear. By means of the tympanometric test AT235 establishes the correct pressure in the ear canal for the subsequent reflex test.

### 6. Testing the next ear

Now select the next ear by pressing "Right" or "Left" (20 or 21) or by pressing the remote switch on the probe. The test procedure is carried out as in the above section 5.

### 7. Print out test result

To print out the test result from the internal printer press "Print" (8).

## Useful Information

- Before starting a test session with tympanometry and reflexes it is possible to select between two pre-programmed reflex setups, namely Test A and Test B. By pressing F6 it is possible to modify the selected Reflex Test A or B, but for the following patient only. Changing the pre-programmed Reflex Test A or B must be done from the Setup, which is reached by pressing F2 while being in the Main Menu.
- Tympanometry and Reflexometry can be performed either on their own or automatically with Reflexometry following a Tympanometry Test. The "Tymp/Reflex" key (10) has two indication lights which indicate whether Tympanometry and Reflexometry is performed on their own or together.
- After an automatic reflex session it is possible to perform manual reflexes simply by pressing "Man" (13). Now, by means of "Frequency Decr/Incr" (14, 15) and "Intensity Decr/Incr" (23, 24), it is possible to alter the reflex parameters. To present the stimuli for the reflex press "Tone Switch" (22). To store the manual reflex press "Store" (16).


## AT235 - Operation Panel Dictionary

| <b>Position:</b> | <b>Symbol:</b> | <b>Function:</b>   |
|------------------|----------------|--|
| 1-6              | F1 – F6        | Function keys related to the functions indicated in the display above the individual function key. |
| 7                | Back           | Pressing "Back" will always take you one level back.   |
| 8                | Print          | Starts printing procedure from internal printer.   |
| 8 + shift        | print screen   | Prints out the present screen / display.   |
| 9                | Pause          | Selects manual start of automatic tympanometry, reflexometry and audiometry tests.                 |
| 9 + shift        | new subject    | Clears the internal memory of the instrument.  |
| 10               | Tymp/Reflex    | Selects between tympanometry, reflexometry and tympanometry + reflexometry.                        |
| 11               | EFT            | Selects Eustachian Tube Function test.   |
| 12               | Audiometry     | Selects audiometry.  |
| 13               | Man            | Select manual tympanometry, reflexometry and audiometry.   |
| 14               | Frequency Decr | Reduces frequency in reflexometry and audiometry.  |
| 15               | Frequency Incr | Increases frequency in reflexometry and audiometry.  |
| 16               | Store          | Stores threshold in manual audiometry and stores manual tympanometry curves and reflexes.          |
| 16 + shift       | no response    | Inserts "no response" sign in the audiogram.   |
| 17               | shift          | Activates the function written in <i>Italic</i> below the buttons 8, 9, 16 and 24.                 |
| 18               | Ipsi           | Selects ipsilateral reflex stimulation.  |
| 19               | Contra         | Selects contralateral reflex stimulation.  |
| 20               | Right          | Selects right ear.   |
| 21               | Left           | Selects left ear.  |
| 22               | Tone Switch    | Presents stimulus in manual reflexometry and audiometry.   |
| 23               | Intensity Decr | Reduces intensity in reflexometry and audiometry.  |
| 24               | Intensity Incr | Increases intensity in reflexometry and audiometry.  |
| 24+shift         | ext. range     | Selects extended range in reflexometry and audiometry.   |

## AT235 - Connection Panel Dictionary

| Position: | Symbol:          | Function:                                  |
|-----------|------------------|--|
| 26        | Left             | Connector for left headphone TDH39.        |
| 27        | Right            | Connector for right headphone TDH39.       |
| 28        | Patient Response | Plug for patient response switch.          |
| 29        | Contra           | Plug for contra telephone.                 |
| 30        | Transducer       | Plug for cable from probe system.          |
| 31        | Air              | Plug for the air hose of the probe system. |
| 32        | Keyboard         | Connection for keyboard                    |
| 33        | Power            | Turns power on/off.                        |
| 34        | USB              | Plug for computer connection.              |
| 35        | Mains            | Plug for mains cable.                      |

### Explanation to symbols which can be found on the instrument:

|   |   |
|---|---|
| I | On (Power: connection to the mains)       |
| O | Off (Power: disconnection from the mains) |
| ~ | Alternating current                       |
|   | Fuse                                      |
|   | Ground                                    |
|   | Dangerous voltage                         |
|   | See explanation in manual                 |
|   | Type B equipment                          |

## Technical Specifications

The technical specifications provided here cover the general aspects of the instrument, while more specific details are found Operation Manual in English.

### Standards:

|            |   |   |
|------------|---|---|
| Audiometer | : | EN 60645-1, ANSI S3.6, Type 4                 |
| Impedance  | : | EN 61027, ANSI S3.39, Type 2                  |
| Safety     | : | EN 60601-1, Class 1, Type B, EMC: EN60601-1-2 |

The CE mark indicates that the instrument is made to fulfil the demands of the Council Directive 93/42/EEC concerning medical devices and that TÜV Product Service, identification no. 0123, has approved the quality system of Interacoustics.

### Impedance:

|                        |                           |
|------------------------|---------------------------|
| Probe Tone Frequency:  | 226 Hz $\pm$ 1%           |
| Probe Tone Intensity : | 85 dB SPL $\pm$ 1.5 dB    |
| Pressure Range :       | +300 to -600 daPa         |
| Accuracy :             | $\pm$ 5% or $\pm$ 10 daPa |
| Compliance Range :     | 0.1-6 ml.                 |
| Compliance Accuracy :  | $\pm$ 5% or 0.1 ml.       |

### Reflex:

|  |  |
|--|--|
| Manual or automatic testing with 18 stimuli to each ear. |  |
| Auto Reflex Detection.                                   |  |
| Multiple Reflex stimuli.                                 |  |
| Stimulus Duration :                                      | 1.0 sec.   |
| Reflex Acceptance :                                      | Adjustable between 2% and 6%, or 0.05 – 0.15 ml change of ear canal volume.    |
| Ipsi Frequencies :                                       | 500, 1000, 2000, 3000, 4000 Hz, wide band, high and low pass.                  |
| Intensity (Max.) :                                       | 100, 105, 110 dBHL.  |
| Contra Frequencies :                                     | 250, 500, 1000, 2000, 3000, 4000, 6000, 8000 Hz, wide band, high and low pass. |
| Intensity (Max.) :                                       | 90, 100, 120, dBHL.  |

### Reflex Decay:

|                      |                             |
|----------------------|-----------------------------|
| Stimulus Duration :  | 10 sec.                     |
| Stimulus Frequency : | } Ipsi and Contra as above. |
| Stimulus Intensity : |                             |

### Audiometer:

Type 4 according to EN 60645-1, ANSI S3.6.

| Frequencies Hz: | Intensities dBHL: |
|-----------------|-------------------|
| 125             | -10 to 90         |
| 250             | -10 to 90         |
| 500             | -10 to 120        |
| 1000            | -10 to 120        |
| 2000            | -10 to 120        |
| 3000            | -10 to 120        |
| 4000            | -10 to 120        |
| 6000            | -10 to 120        |
| 8000            | -10 to 110        |

Auto Threshold Determination (Modified Hughson Westlake).

# AT235 Instruction for Use - English

Date: 2001-06-15

Page 6/7

## Tests:

|                                  |   |  |
|----------------------------------|---|--|
| <i>Tympanometry</i>              | : | Automatic  |
| <i>Reflexometry incl. decay:</i> | : | Automatic or Manual.   |
| <i>Audiometry</i>                | : | Manual or Auto Threshold (Patient controlled Hughson Westlake test after ISO 8253) |
| <i>Eustachian Tube</i>           | : | Semi Automatic function test.<br>Williams Test (non perforated eardrums)           |

## Calibration:

Impedance : EN 61027, ANSI S3.39

Audiometer : ISO 389-1, ANSI S3.6

Calibration is performed via the instrument's front panel and is stored in a permanent memory.

## Computer Communication:

Built-in USB input/output computer interface. One mode allows the computer to both monitor and control the AT235. Control actions will be indicated on the operation panel and the display.

## Power supply:

100, 110, 120 V AC  $\pm 10\%$ . Fuses: T 1A

220, 230, 240 V AC  $\pm 10\%$ . Fuses: T 0.5A

50 - 60 Hz.

100 VA maximum

## Dimensions:

W x D x H: 36 x 26 x 10 cm, 14 x 10 x 4 inches. Weight: 2.8 kg / 6 lbs.

## Included Parts

- AT235-ATP Standard Probe System
- TDH39 Single Contralateral Earphone
- UPS400 External Switch Mode Power Supply
- Power Cable (110/220V, please specify)
- BET50 Box of 65 Assorted Eartips
- TPR35 3 Rolls of Recording Paper
- PCR-AT235 Dust Cover
- Operation/Multilingual CE manuals

### Concerning repair

1. Interacoustics is only considered to be responsible for the validity of the CE marking, effects on safety, reliability and performance of the equipment if:

assembly operations, extensions, readjustments, modifications or repairs are carried out by authorised persons,

a 1 year service interval is maintained


the electrical installation of the relevant room complies with the appropriate requirements, and

the equipment is used by authorised personnel in accordance with the documentation supplied by Interacoustics.

2. It is important that the customer (agent) fills out the RETURN REPORT every time a problem arises and sends it to Interacoustics, Drejervaenget 8, DK-5610 Assens, Denmark. This should also be done every time an instrument is returned to Interacoustics. (This of course also applies in the unthinkable worst case of death or serious deterioration to patient or user).

3. When instrument fuses need renewal, the correct type as stated on the instrument shall be used.

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